A Parser for mzXML, mzData and mzML files

Bernd Fischer, Steffen Neumann

December 19, 2010

Contents

\$analyzer
[1] "FTMS"

```
A short example sequence to read data from a mass spectrometer. First open the file.
> library(mzR)
> library(msdata)
> mzxml <- system.file("threonine/threonine_i2_e35_pH_tree.mzXML",
                        package = "msdata")
> ramp <- rampOpenFile(mzxml)</pre>
We can obtain different kind of header information.
> ramp$getRunInfo()
$scanCount
[1] 55
$lowMZ
[1] 1.58101e-322
$highMZ
[1] 2.371515e-322
$startMZ
[1] 1.513193e-316
$endMZ
[1] 6.928674e-310
$dStartTime
[1] 0.3485
$dEndTime
[1] 390.027
> ramp$getInstrumentInfo()
$manufacturer
[1] "Thermo Scientific"
$model
[1] "LTQ Orbitrap"
$ionisation
[1] "ESI"
```

\$detector

[1] "unknown"

- > SH <- ramp\$getAllScanHeaderInfo()</pre>
- > head(SH)

	seqNum acc	quisitionNum msLevel	peaksCou	int tot	IonCurrent	retenti	onTime
1	1	1 1	. 6	884	341427000		0.3485
2	2	2 2	2 4	32	160473000		5.8561
3	3	3 2	2 3	340	58862000	1	2.5000
4	4	4 3	3 2	273	30770400	1	9.7568
5	5	5 3	3 2	238	5291800	2	26.6056
6	6	6 3	3 1	.40	14294000	3	34.0324
	basePeakM2	Z basePeakIntensity	collision	Energy	ionisation	nEnergy	lowMZ
1	120.0660	211860000		0		0	50.3254
2	120.0660	139169000		0		0	50.4459
3	102.0560	27036600		35		0	50.0658
4	102.0550	26736000		0		0	50.2843
5	56.0497	2188950		35		0	50.1967
6	74.0605	12518700		0		0	50.2254
highMZ precursorScanNum precursorMZ precursorCharge precursorIntensity							
1	298.6730	0	0.0000		0		0
2	134.1380	0	120.0661		1		210140000
3	134.1430	0	120.0661		1		210140000
4	114.6510	0	102.0600		0		1441880
5	114.9470	0	102.0600		0		1441880
6	84.8802	0	74.0600		0		1556090
${\tt mergedScan}\ {\tt mergedResultScanNum}\ {\tt mergedResultStartScanNum}$							
1	()	0		(0	
2	()	0		(0	
3	()	0		(0	
4	()	0			0	
5	()	0		(0	
6	C		0		(0	
	mergedResultEndScanNum						
1		0					
2		0					
3		0					
4		0					
5		0					

Read a single spectrum from the file.

0


- > PL <-ramp\$getPeakList(10)</pre>
- > PL\$peaksCount

[1] 317

6

> head(PL\$peaks)

[,1] [,2]
[1,] 50.08176 6984.858
[2,] 50.62267 7719.419
[3,] 50.70530 7185.290
[4,] 50.73298 7509.140
[5,] 50.83848 9366.624
[6,] 50.88303 8012.808



You should close the file when not needed any more. This will release the memory of cached content.

> ramp\$close()